

Cycle sequencing protocol

For each sample, the target gene will be sequenced in both forward and reverse directions. Therefore, we need to prepare **two separate cycle sequencing reactions per sample**.

Osea, dos placas.

1. Make BigDye mix.

	X1	X72	X81
Nuclease-free H ₂ O	6.25µl	450 µl	506.25 µl
BigDye Buffer	1.75µl	126 µl	141.75 µl
10 µM Forward OR Reverse primer	0.50 µl	36 µl	40.5 µl
Primer BigDye	0.50µl	36 µl	40.5 µl
DNA Sample	1.00µl		
TOTAL	10 µl	648 µl	729 µl

** BigDye is light sensitive, so keep it in freezer until you need to add.

2. In a clean 96-well plate, add to each well:

9 µl BigDye mix + 1 µl cleaned PCR product (from SPRI bead clean-up)

3. Run in thermocycler:

95°C for 0:30 sec	35x
50°C for 0:30 sec	
60°C for 4:00 min	
10°C	forever

** Upon completion, BigDye plates must be stored in the refrigerator (4°C), **not in** the freezer.

Protocol

1. Prepare two (**2**) **BigDye mix** tubes, one for the forward reaction and one for the reverse reaction. Both of X81 volume.
2. Agregas 9 µl de BigDye mix a cada pocillo (72) y 1 µl del clean PCR product.
3. Corres en el termociclador con el programa guardado como BigDye.