

## GENEIOUS SEQUENCING ALIGNMENT, CLEANING, AND TREE BUILDING

1. Import sequences
  - a. Open geneious and import sequences to a "raw sequences" folder.
  - b. Check the percentage of untrimmed bases from each sequence (%HQ) > 30% is good, >50% is optimal.
2. Alignment *de novo* to generate contigs (set of DNA segments or sequences that overlap in a way that provides a contiguous representation of a genomic region).
  - a. Select pairs of sequences to align
  - b. Click on "align/assemble"
  - c. Click on "de novo assemble"
    - i. Options Data box -> Assemble by "1<sup>st</sup>" part of name, separated by "\_ (Underscore)", or choose as appropriate. Keep all other options as default.
3. Save aligned sequences into folder "aligned sequences"
4. Clean aligned sequences by:
  - a. Delete regions with pink bars if "small". If pink regions are too large check and assess whether delete a large section is feasible.
    - i. Select "allow editing" to delete pink sections
    - ii. If pink trim regions are not annotated already, select "annotate & predict" next to "allow editing"
    - iii. Select "trim ends" keep default values -> then delete pink trim regions
  - b. Check gaps and ambiguities
    - i. On the right icon tabs, select the monitor icon -> highlighting -> drop down select gaps/ambiguities -> go <> "in any sequence"
    - ii. On ambiguities or gaps highlighted sections, check line graph with highest peak (Red-Adenine, Blue-Cytosine, Green-Thymine, Yellow-Guanine) to make informed decisions about base changes/filling gaps.
    - iii. Save. Check "NO" when asked *"Do you wish to apply changes to the original (unaligned) sequences?"*  
*If you choose no, this alignment will permanently discard references to the unaligned sequences and you will no longer be able to apply changes to them"*
    - iv. Copy and paste clean contigs into a new folder named "clean contigs".
5. Perform a multiple alignment with all clean contigs
  - a. Click on align/assemble -> multiple align -> consensus align -> create alignment of consensus sequences only
  - b. Copy clean contigs alignment into a new folder "alignment genebank"
6. Perform multiple alignment with genebank sequences
  - a. Download sequences from genebank (find relevant sequences from the literature) into the "alignment genebank" folder
  - b. Select genebank sequences and clean contig alignment
  - c. Click on align/assemble -> consensus align -> create alignment of all sequences-> check box "automatically determine direction"

- d. Trim aligned sequences into smallest fragment
- 7. Copy genebank alignment into new folder "tree"